BEC SARS-CoV-2 RT-LAMP kit for human samples

BEC SARS-CoV-2 RT-LAMP kit for human samples

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English

► **Intended Use**

The BEC SARS-CoV-2 RT-LAMP kit for human samples is a test for rapid detection of SARS-CoV-2 combining a technology of rapid extraction of RNA from nasal or nasopharyngeal samples and a RT-LAMP molecular test targeting 2 genes of SARS-CoV-2.

Positive results are indicative of the presence of SARS-CoV-2 RNA; clinical correlation with patient history and other diagnostic information is necessary to determine patient infection status.

Negative results do not preclude SARS-CoV-2 infection and should not be used as the sole basis for patient management decisions. Negative results must be combined with clinical observations, patient history and epidemiological information.

The analysis of this test must be carried out with the BEC Reader device.

► **Content of the kit**

Number of tests in a kit : 10.

Kit storage temperature : +4°C.

Kit shipping temperature: room temperature (RT) or +4°C.
The kit is stable at RT for 2 weeks.

Kit shelf life : expiry date indicated on the outside of the box.

Kit content :

- 5 individual pouches containing 2 tests each.
- A cardboard rack, which makes it easier to handle the tubes, see Figure 2.
- A user guide.

Content of a pouch :

Description	Reference	Color	Quantity per pouch	Volume
BEC SARS-CoV-2 Assay Tubes for human samples	P04301	-	1 strip of 6 microtubes + 2 strips of caps	-
BEC SARS-CoV-2 Lysing Buffer (green)	P13300	Green	2	1.050 mL
BEC SARS-CoV-2 Neutralization Buffer (blue)	P20300	Blue	2	1 mL
BEC SARS-CoV-2 Reaction Buffer (orange)	P21300	Orange	2	150 μ L
BEC SARS-CoV-2 Transfer Tubes (violet)	P30300	Violet	2	-
sterile 3mL syringe	P31300	-	2	-
sterile needle for 3mL syringe	P32300	-	2	-
5 micron filter for 3mL syringe	P33300	-	2	-

► **Materials required but not provided**

- The BEC field reader pack including :
 - The BEC Reader device (see Picture 2)
 - A 25 µL fixed-volume pipette (orange)
 - A 100 µL fixed-volume pipette (violet)
 - A box of disposable filter-tips (recommended CLEARLINE cat. number 713116)
 - Biological waste bin
- Sampling material : Swabs for nasal or nasopharyngeal sample
Examples of swabs:
 - iClean® Specimen Collection Flocked Swab
 - Dry Swabs 155C, COPAN Diagnostics
 - Dry Swabs 159C, COPAN Diagnostics
- Personal protective equipment for the user according to the local regulation.

► **Warnings and precautions**

Please read the user manual carefully before using the product.

- Before any use, please check that the product and its contents are not deteriorated and that all the contents are in the kit.
- Do not use a test from a pouch that has been opened for more than 24 hours, regardless of the storage temperature.
- Handle the samples as if they were infectious and/or dangerous, in accordance with the security procedures applicable in the laboratory.
- Wear adequate personal protective equipment (laboratory coat, disposable powder-free gloves,

mask).

- The operator must handle the syringe with the needle with care. To avoid any risk of injury, the needle protective cap must be put back on it between each preparation step. At the end of handling, the needle with its protective cap must be discarded in accordance with local safety rules. Avoid microbiological and amplicon contaminations (by DNase/ RNase) of the sample and of the kit contents.
- Always use pipettes with disposable DNase/RNase-free filter-tips
- The user must check the accuracy of the pipettes regularly
- Always wear protective powder-free gloves when handling the contents of the kit.
- Do not open the reaction tubes after amplification to avoid any contamination with amplicons.
- Do not use the kit content after their expiry date.
- Dispose of the samples, contents in the pouch and waste material resulting from the test in accordance with local safety regulations.
- This kit is intended for use with the reader approved for this purpose. Any directly interpretable result is dependent on the reader settings. If another reading instrument is used, it is the responsibility of the operator to ensure the correct calibration of that instrument. Bertin Technologies shall bear no responsibility for results obtained with any other reference system.
- Do not leave the microtubes in the reader once the reaction is completed.
- It is recommended to use a flat, clean and clear working surface in a covered and closed area at a temperature between 15°C and 25°C.
- Before opening the preparation tubes, make sure that all the liquid or lyophilized pellet contained in

the tube is at the bottom of the tube. If needed, tap the tube several times downwards to lower any drops.

Never reopen the microtubes once they have been inserted in the BEC Reader device, nor after the end of the reaction. Re-opening under these conditions would result in environmental contamination. This contamination may cause false positive results for subsequent tests.

► **General information**

The BEC SARS-CoV-2 RT-LAMP kit for human samples is a « Point-of-Care » molecular diagnostic kit allowing the detection of the SARS-CoV-2 virus in nasal or nasopharyngeal samples.

The SARS-CoV-2 is a coronavirus responsible for the emergence of a new respiratory infectious disease called Covid-19 that has caused a global pandemic. This disease, which appeared in China at the end of 2019, is transmitted mainly by air, with a high reproduction rate of between 2 and 3 in the absence of control and prevention measures. The severity of the clinical signs of the disease varies greatly from one individual to another, with a preponderance of severe forms in elderly individuals or those suffering from other diseases.

The scientific data available to date make it possible to estimate that 30% to 60% of infected subjects are asymptomatic but can transmit the disease. In this context, it is essential to be able to carry out diagnostic tests on a massive scale in order to limit the spread of the virus.

The BEC-SARS-CoV-2 kit has been developed in order to allow diagnostic tests to be carried out by non-specialized personnel in remote health establishments or on dedicated test sites.

Source : Covid-19 Disease Sheets – Institut Pasteur – May 2020 (<https://www.pasteur.fr/en/medical-center/disease-sheets/covid-19-disease-novel-coronavirus>)

► **Product description**

The BEC SARS-CoV-2 RT-LAMP kit for human samples is an ensemble of reagents allowing rapid extraction of viral RNA and real-time RT-LAMP amplification of 2 sequences of the SARS-CoV-2 virus genome. It contains an endogenous internal control to ensure the correct running of the extraction and amplification.

RT-LAMP is a method of specific isothermal amplification of a nucleotide sequence with incorporation of a fluorescent nucleic acid stain during amplification (Li & Macdonald, 2015; Nagamine et al., 2002; Notomi et al., 2000; Yan et al., 2014). It amplifies a target sequence with two to three sets of primers at a constant temperature of 65°C. Compared to traditional detection methods, RT-LAMP is more specific and faster (Aebischer et al., 2014). It produces more amplicons and is less sensitive to inhibitors that are naturally present in biological matrices (Chander et al., 2014; Francois et al., 2011; Nie et al., 2012).

The BEC SARS-CoV-2 RT-LAMP kit for human samples allows 2 samples to be tested in parallel.

An analysis with the BEC SARS-CoV-2 RT-LAMP kit for human samples consists of carrying out the following steps :

- For the 2 samples in a sequential manner :
 - o Rapid extraction of viral RNA with chemical lysis.
 - o Neutralization of the amplification inhibitors present in the lysate by precipitation with chaotropic agents.
 - o Filtration of the lysate to remove aggregates of cell debris and reagents precipitated during neutralization.
- For the 2 samples in parallel: RT-LAMP amplification with the reaction buffer and the primers distributed in the microtubes (assay tubes), see Figure 1.


Figure 1: A strip of 6 microtubes allowing the amplification of 3 sequences of the RdRP, N and Rnase P genes (CI : internal control) for 2 samples

The test allows the amplification of 3 targeted sequences in the following genes:

- RdRP Gene : target specific to betacoronavirus (β CoV) human pathogens belonging to B and C groups (SARS-CoV-2, SARS-CoV and MERS-CoV) - positions 2 and 5 on the strip of microtubes (see Figure 1)
- N Gene : target specific to SARS-CoV-2 - positions 3 and 6 on the strip of microtubes (see Figure 1)
- RNase P, human gene used as internal control for the extraction and amplification processes - positions 1 and 4 on the strip of microtubes (see Figure 1).

These amplifications are performed in parallel and the amplification result is obtained after a period of between 5 and 25 minutes depending on the viral load of the sample being analyzed.

The reaction buffer used to perform RT-LAMP amplification contains a 6-FAM fluorescent dye. If the virus is present in the sample, SARS-CoV-2 cDNA is amplified and 6-FAM fluorescent nucleic acid stains are incorporated into the generated nucleotide sequences. The BEC Reader device measures the fluorescent signal emitted on the FAM channel by the sample every ten seconds for 25 minutes. An algorithm analyses the slope and the amplitude of fluorescent signal to determine if there is a significant amplification.

The BEC Reader device indicates in real time for each tube whether the amplification is positive (detection of the RNA target) or negative (no detection of the RNA target).

► Operating Protocol

▷ Preparation of the samples

The test is compatible with nasal and nasopharyngeal samples collected with a swab.

For nasopharyngeal sample, the following protocol is recommended:

- Tilt patient's head back 70 degrees.
- Gently and slowly, insert the swab through the nostril parallel to the palate (not upwards) until resistance is encountered.
- Gently rub and roll the swab.
- Leave swab in place for several seconds to absorb secretions.
- Slowly remove swab while rotating it.
- Repeat the protocol in the second nostril with the same swab

For nasal sample, the following protocol is recommended:

- Insert the swab into the nostril (1-2 cm)
- Gently rotate 5 times against the walls of the nostril.
- Repeat the protocol in the other nostril with the same swab

The sample must be stored and/or transported in a dry form at +4°C and analyzed within 24 hours after the sampling

CAUTION! The test shows a significant decrease in sensitivity when used with samples previously placed in a liquid medium between collection and testing.

▷ Reagent preparation

- Open the aluminum pouch and check that the buffers or lyophilized pellets contained in the various tubes are well located at the bottom of the tubes and that no tubes are open. If this is the case, do not use the pouch.
If necessary, tap the tubes on the benchtop to make sure that drops of buffer or lyophilized pellets do not remain in the cap.
- Leave the reagents at room temperature for 5 minutes before use. If the solution contained in the green tube is crystallized, warm the tube in your hands for a few seconds.
- Place the different tubes on the cardboard rack provided with the kit, respecting the color codes, see Figure 2. **Do not open the tubes at this stage.**


Figure 2: Rack that is provided for easy handling of the tubes.

▶ **BEC Reader device preparation**


- Switch on the BEC Reader device by pushing  button in front of the device during 3 seconds. The screen turns on and the reader performs a self-test. The cover opens.


Figure 3: BEC Reader device with location of ON/OFF switch button.


- Press on the screen and then press on TEST symbol


- Select the test "BEC_Human Assay Vx.x"


- Complete the RunID : this is the identifier of the run result file.
- Press "Auto 1-8" to automatically assign a name to the tubes. Go to the next screen by pressing ✓

- From this stage, the reader goes into preheating until reaching the set temperature of 65 °C.

- When the setpoint is reached, the reader informs the user.

▷ **Test execution**

CAUTION : perform the steps 1 through 5 with the 1st sample, and then repeat these steps with the 2nd sample. If only one sample is available, perform steps 1 through 6 and leave the wells 4 through 6 empty.

1) Step 1 : lysis and RNA extraction

- a) Open the green tube (BEC SARS-CoV-2 Lysing Buffer (green)) in position 1 on the rack (see Figure 2).
- b) Immerse the swab in the solution and turn and stir the swab in the solution.
- c) Squeeze the swab against the inner side of the tube and discard the swab.
- d) Close the green tube, shake it vigorously during 10 seconds and place it in position 2 on the rack.

DO NOT leave the sample to stand before proceeding to Step 2 to avoid degradation of the mix that would interfere with the test.

2) Step 2 : neutralization

- a) Open the blue tube (BEC SARS-CoV-2 Neutralization Buffer (blue)) and the green tube in position 2 on the rack.
- b) Withdraw the entirety of the liquid from the blue tube with the syringe fitted with the needle and dispense it in the green tube.
DO NOT pour directly without using the needle.
- c) Lay the syringe aside.
- d) Close the blue tube and the green tube. Shake the

green tube vigorously during 10 seconds and place it on the rack in position 3.

3) Step 3 : filtration

- a) Tap the green tube on the benchtop to lower the droplets that are in the cap. Open the green tube and the violet tube (BEC SARS-CoV-2 Transfer Tubes (violet)) in position 3.
- b) Draw as much liquid as possible from the green tube using the syringe fitted with the needle. Be careful to draw liquid from the green tube and not foam.
- c) Take out the needle from the syringe and replace it with the filter. Discard the needle in an appropriate waste bin.
- d) Insert the filter into the violet tube. Gently push the syringe plunger until it stops. By gently pushing the plunger, the flow of the filtrate is facilitated.
- e) Close the green tube and the violet tube. Place the violet tube in position 4 on the rack. Discard the green tube in an appropriate waste bin.

4) Step 4 : transfer

- a) Open the violet tube and the orange tube (BEC SARS-CoV-2 Reaction Buffer (orange)) in position 4 on the rack.
- b) Aspirate a 100 μL volume from the violet tube with the violet fixed-volume pipette fitted with a filter-tip and transfer it to the orange tube.
- c) Close the violet tube and the orange tube. Shake vigorously the orange tube during 10 seconds and place it in position 5 on the rack.

CAUTION : Orange tube should be used in step 5 and the incubation (step 6) should be started within 30 minutes after preparation of orange tube.

5) Step 5 : dispensing in the microtubes

a) For sample #1

- Open the orange tube in position 5.
- Open the microtubes of the strip and discard the strip of caps.
- Use the orange pipette fitted with a filter-tip to aspirate a 25 μ L volume from the orange tube and dispense this volume in tubes # 1, 2 and 3 (see Figure 1).

Use a new filter tip for each microtube

- Close the orange tube.
- Discard the green, blue, violet and orange tubes

CAUTION : Repeat steps 1 through 5 with the second sample

b) For sample #2

- Open the orange tube in position 5.
- Use the orange pipette fitted with a filter-tip to aspirate a 25 μ L volume from the orange tube and dispense this volume in tubes # 4, 5 and 6 (see Figure 1).

Use a new filter tip for each microtube

- Close the orange tube.
- Use the additional strip of caps to close the microtubes of the strip.
- Discard the green, blue, violet and orange tubes.

c) When the strip of microtubes is filled and closed, proceed immediately to step 6.

6) Step 6 : Incubation



CAUTION : Make sure the BEC reader has reached the temperature set point (65°C).

- Shake the microtubes until the lyophilized pellets are completely dissolved. Then tap the tubes on the benchtop to lower all the droplets to the bottom of the tube.
- Insert the microtube strip into the BEC Reader device, placing the black end tab on the left side of the instrument (see Figure 4). Leave the wells 7 and 8 empty.

Figure 4: How to insert microtube strip in the BEC Reader device. Red arrows showing position of the black end tab. Blue arrows showing empty well positions 7 and 8.

- c) Immediately close the lid. The BEC Reader device locks and starts the program automatically.
- d) A progress bar allows you to follow the progress of the test. The result for each well is displayed in real time on the screen. At the end of the amplification, the message "test completed" appears and all the results are displayed.

- e) The lid opens automatically at the end of the amplification.
- f) Discard immediately the microtubes **without opening the microtubes**.

► Interpretation of results

► Pictograms in positions 1 through 3 and 4 through 6

For each of the two samples (see Figure 1 for the position), the analysis of the results shall be performed pursuant to the following table :

Internal Control	RdRP gene	N gene	Interpretation of results
			Detection of viral RNA specific to the betacoronavirus groups B and C (SARS-CoV-2, SARS-CoV-1 and MERS-CoV) and to the SARS-CoV-2 virus Sample is SARS-CoV-2 positive
			Detection of viral RNA specific to betacoronavirus groups B and C (SARS-CoV-2, SARS-CoV-1 and MERS-CoV) No detection of viral RNA specific to the SARS-CoV-2 Sample is presumptive positive for SARS-CoV-2 Repeat the test with a new sample.*
			Detection of viral RNA specific to SARS-CoV-2 virus No detection of viral RNA specific to betacoronavirus groups B and C (SARS-CoV-2, SARS-CoV-1 and MERS-CoV) Sample is SARS-CoV-2 positive

Internal Control	RdRP gene	N gene	Interpretation of results
			No detection of viral RNA specific to betacoronavirus groups B and C (SARS-CoV-2, SARS-CoV-1 and MERS-CoV) and to SARS-CoV-2 virus Sample is SARS-CoV-2 negative
			Detection of viral RNA specific to betacoronavirus groups B and C (SARS-CoV-2, SARS-CoV-1 and MERS-CoV) and to SARS-CoV-2 virus No detection of RNase P human gene Sample is SARS-CoV-2 positive
			Detection viral RNA specific to betacoronavirus groups B and C (SARS-CoV-2, SARS-CoV-1 and MERS-CoV) No detection of viral RNA specific to SARS-CoV-2 No detection of RNase P human gene Sample is presumptive positive for SARS-CoV-2 Repeat the test with a new sample.*
			Detection of viral RNA specific to SARS-CoV-2 virus No detection of viral RNA specific to betacoronavirus groups B and C (SARS-CoV-2, SARS-CoV-1 and MERS-CoV) No detection of RNase P human gene Sample is SARS-CoV-2 positive
			RT-LAMP inhibition or reagent failure. Repeat the test with a new sample.**

*: If the results of the initial test and repeat test are identical, use an alternative method for SARS-CoV-2 molecular detection such as RT-qPCR with a new sample.

** : If the results of the initial test and repeat test are identical, please contact technical support.

Note : The « unclear »  pictogram may appear on one or several of the positions 1 through 6. In this case, the test is not valid as the amplification did not proceed properly. The test must be repeated with a new sample. In case the « unclear » pictogram appears on two successive runs, contact Technical Support.

▷ **Pictograms in positions 7 and 8**

Positions 7 and 8 always show the « unclear »  pictogram.

▶ **Optional: exporting results**

The BEC reader device automatically saves the test results into its memory up to 50 tests.

The test results can be exported as files in .CSV and .pdf format. Both files can be opened on a PC.

Pdf files contain run identification information (test name, RunID, date and time) and the results that are displayed on the screen.

To access detailed results (including fluorescence data), .CSV files must be opened with the 'T8 ISO Desktop' software.


- Connect the provided USB stick in the plug in front of the device

- Press the "Results" icon on the Home Screen. The list of test results is displayed.

- Select the test results and press the Export Results icon  to export results to the USB key.

Be sure the export is completed before removing the USB stick.

► Limitations

- Performance was evaluated with nasal and nasopharyngeal swab samples only, following the procedures described below. Failure to follow these procedures may result in assay failure or cause erroneous results.
- The test is a qualitative test and does not provide the quantitative value of detected organism present.
- This test cannot rule out diseases caused by other bacterial or viral pathogens.
- Interpretation of results should account for the possibility of false negative and false positive results (see table below).

Table: Potential causes of false negative and false positive results

False negative results	False positive results
Improper collection, handling and/or storage of samples	Improper handling of samples
Samples with a viral load below the limit of detection	Contamination of workspaces
Improper sample preparation before analysis	Opening of reaction tubes post-amplification
Virus mutations in the regions targeted by the test	Deviation from handbook protocol
The presence of RT-LAMP inhibitors or interfering substances	
Deviation from handbook protocol	

► **Quality Control**

In accordance with Bertin Technologies Quality Management System (ISO 9001:2015 and ISO 13485:2016), critical components of the kit are tested against predetermined specifications to ensure consistent product quality. Quality control testing is performed using standard templates with results compared to previous lots.

► **Performance characteristics**

▷ **Clinical performance (nasal swab)**

Clinical performance characteristics of BEC SARS-CoV-2 RT-LAMP kit for human samples was evaluated on nasal swab in a prospective study in Guadeloupe. Symptomatic and asymptomatic patients were included for this study.

For each patient, one nasal swab and one nasopharyngeal swab were performed. The nasal swab was tested with the BEC SARS-CoV-2 RT-LAMP kit for human samples according to product instructions and the nasopharyngeal swab was eluted in viral transport media (VTM) and tested using the TaqPath™ COVID-19 CE-IVD RT-PCR kit from ThermoFisher Scientific (3 targeted viral regions of SARS-CoV-2 : N, ORF1b and S).

The diagnostic sensitivity and specificity were determined using a panel of 71 individual nasal swabs by comparison with RT-qPCR results obtained on nasopharyngeal swabs. Of the 71 patients, 23 patients presented a negative RT-qPCR test and 48 patients presented a positive RT-qPCR test.

	BEC SARS-CoV-2 RT-LAMP kit
% Detection of samples with $C_T \leq 20^*$ (10 qRT-PCR positive samples)	100%
% Detection of samples with $C_T \leq 25^*$ (21 qRT-PCR positive samples)	91,9%
% Detection of samples with $C_T \leq 30^*$ (34 qRT-PCR positive samples)	77,9%
% Detection of samples with $C_T \leq 33^*$ (43 qRT-PCR positive samples)	64,3%
% Detection of samples with $C_T \leq 40^*$ (48 qRT-PCR positive samples)	58,3%
Overall Diagnostic Sensitivity (Se) (Concordance with RT-qPCR to detect positives)	100% (samples with $C_T \leq 20$) 91,9% (samples with $C_T \leq 25$) 77,9% (samples with $C_T \leq 30$) 64,3% (samples with $C_T \leq 33$) 58,3% (samples with $C_T \leq 40$)
Overall Diagnostic Specificity (Sp) (Concordance with RT-qPCR to detect negatives)	100 %

* C_T values correspond to equivalent C_T values obtained with IP2 and IP4 RT-qPCR reference method.

▷ **Clinical performance (nasopharyngeal swab)**

Clinical performance characteristics of BEC SARS-CoV-2 RT-LAMP kit for human samples was evaluated on nasopharyngeal swab in a prospective study in the Parisian hospital La Pitié Salpêtrière. The study was carried out on patients hospitalized in the infectious disease department due to the covid-19. All patients enrolled in the study tested positive for covid-19 (RT-qPCR) upon admission.

For each patient, two nasopharyngeal swabs were performed. One swab was tested with the BEC SARS-CoV-2 RT-LAMP kit for human samples according to product instructions and the other one was eluted in viral transport media (VTM) and tested using the TaqPath™ COVID-19 CE-IVD RT-PCR kit from ThermoFisher Scientific (3 targeted viral regions of SARS-CoV-2 : N, ORF1b and S). All samples were analyzed within 3 hours after the sampling.

The diagnostic sensitivity and specificity were determined using a panel of 58 swabs by comparison with RT-qPCR results. Of the 58 patients, 22 patients presented a negative RT-qPCR test and 36 patients presented a positive RT-qPCR test.

	BEC SARS-CoV-2 RT-LAMP kit
% Detection of samples with $C_T \leq 20^*$ (0 qRT-PCR positive sample)	NA
% Detection of samples with $C_T \leq 25^*$ (2 qRT-PCR positive samples)	100%
% Detection of samples with $C_T \leq 30^*$ (13 qRT-PCR positive samples)	81,6%
% Detection of samples with $C_T \leq 33^*$ (26 qRT-PCR positive samples)	62,6%
% Detection of samples with $C_T \leq 40^*$ (36 qRT-PCR positive samples)	52,4%
Overall Diagnostic Sensitivity (Se) (Concordance with RT-qPCR to detect positives)	100% (samples with $C_T \leq 25$) 81,6% (samples with $C_T \leq 30$) 62,6% (samples with $C_T \leq 33$) 52,4% (samples with $C_T \leq 40$)
Overall Diagnostic Specificity (Sp) (Concordance with RT-qPCR to detect negatives)	95 % (1 false positive)

* C_T values correspond to equivalent C_T values obtained with IP2 and IP4 RT-qPCR reference method.

▷ Analytical performance

Analytical performance of the BEC SARS-CoV-2 RT-LAMP kit has been evaluated by the Laboratory for Urgent Response to Biological Threats from the Environment and Infectious Risks Research and Expertise Unit of Pasteur Institute.

■ Analytical sensitivity

To determine the analytical sensitivity of the BEC RT-LAMP reaction, serial 10-fold dilutions of genomic SARS-CoV-2 were tested in parallel with the reference real-time RT-PCR developed by the National Reference Centre (NRC) for Respiratory Viruses in France (WHO Coronavirus disease COVID-19 technical guidance: Laboratory testing for 2019-nCoV in humans, available from

https://www.who.int/docs/default-source/coronaviruse/whoinhouseassays.pdf?sfvrsn=de3a76aa_2) and the BEC RT-LAMP reaction.

Each dilution of template RNA was tested 8 times.

Table: Analytical sensitivity (limit-of-detection) of the BEC RT-LAMP reaction as compared to the gold-standard RT-PCR assay

Assay	Target gene	Primer set ID	Quantified genomic SARS-CoV-2 RNA					
			10 genome copies			100 genome copies		
			Mean Cp	SD Cp	Detected replicates (out of 8)	Mean Cp	SD Cp	Detected replicates (out of 8)
Real-time RT-PCR								
<i>RdRp</i>		IP2 (FAM)	37.02	0.56	8	34.65	0.48	8
		IP4 (HEX)	38.18	0.86	7	35.05	0.36	8
Real-time RT-LAMP								
<i>N</i>		E	11.48	2.66	8	7.95	0.38	8
		<i>RdRp</i>	T	11.21	2.47	7	7.59	0.71

The BEC RT-LAMP reaction demonstrated a limit of detection of 10 genome copies per reaction.

Clinical sensitivity of the BEC RT-LAMP reaction was then evaluated on 62 extracts of respiratory samples from SARS-CoV-2 positive patients.

		RT-qPCR reference method results	
		POS	NEG
BEC RT-LAMP reaction results	POS	59	0
	NEG	3	53
Total		62	53

Clinical sensitivity of the BEC RT-LAMP reaction is 95%.

- Analytical specificity

The BEC SARS-CoV-2 RT-LAMP kit allows the amplification of 2 targeted sequences in the following genes:

- RdRP Gene : target specific to betacoronavirus (β CoV) human pathogens belonging to B and C groups (SARS-CoV-2, SARS-CoV and MERS-CoV)
- N Gene : target specific to SARS-CoV-2

Analytical specificity of the BEC SARS-CoV-2 RT-LAMP kit was validated by *in silico* analysis of the primers in comparison with all SARS-CoV-2 and β CoV B and C groups sequences available in the worldwide databases (April 2020).

A complementary *in silico* analysis of the primers was performed in January 2021 to evaluate primers capability to detect both English and South African variant strains of SARS-CoV-2. The results of this analysis predicted that BEC SARS-CoV-2 RT-LAMP kit is able to amplify both variant strains. Comparative sequences analysis identified only one mismatch located on the first nucleotide of 5' end of one in six primers (F3) for RdRP detection. Due to its location on 5' end, this mismatch is not predicted to alter performance amplification.

Analytical specificity was evaluated experimentally using a panel of extracts of RNA from human pathogens viruses:

- Coronavirus OC43 (1 extracted strain)
- Coronavirus 229E (1 extracted strain)
- Influenza A/H1N1 (1 extracted strain)
- Influenza A/H3N2 (1 extracted strain)
- Influenza B (1 extracted strain)

No cross reactivity of the BEC SARS-CoV-2 RT-LAMP kit with the selected microorganisms was observed.

Analytical specificity was also evaluated with 45 extracts of respiratory samples from SARS-CoV-2 negative patients.

No false positive was observed.

► Technical Support

Before you contact Bertin Bioreagent Technical Support, please gather the following information :

- Product Name
- Batch Number
- Serial Number of the instrument
- Error messages (if any)

Technical Support contact details :

- +33 (0)139 306 036
- tech@bertin-bioreagent.com

► Bibliography

Aebischer, A., Wernike, K., Hoffmann, B., & Beer, M. (2014). Rapid genome detection of Schmallenberg virus and bovine viral diarrhoea virus by use of isothermal amplification methods and high-speed real-time reverse transcriptase PCR. *Journal of Clinical Microbiology*, 52(6), 1883–1892. <https://doi.org/10.1128/JCM.00167-14>

Chander, Y., Koelbl, J., Puckett, J., Moser, M. J., Klingele, A. J., Liles, M. R., Carrias, A., Mead, D. A., & Schoenfeld, T. W. (2014). A novel thermostable polymerase for RNA and DNA loop-mediated isothermal amplification (LAMP). *Frontiers in Microbiology*, 5(AUG). <https://doi.org/10.3389/fmicb.2014.00395>

Francois, P., Tangomo, M., Hibbs, J., Bonetti, E. J., Boehme, C. C., Notomi, T., Perkins, M. D., & Schrenzel, J. (2011).

Robustness of a loop-mediated isothermal amplification reaction for diagnostic applications. *FEMS Immunology and Medical Microbiology*, 62(1), 41–48. <https://doi.org/10.1111/j.1574-695X.2011.00785.x>

Li, J., & Macdonald, J. (2015). Advances in isothermal amplification: Novel strategies inspired by biological processes. In *Biosensors and Bioelectronics* (Vol. 64, pp. 196–211). Elsevier Ltd. <https://doi.org/10.1016/j.bios.2014.08.069>

Nagamine, K., Hase, T., & Notomi, T. (2002). Accelerated reaction by loop-mediated isothermal amplification using loop primers. *Molecular and Cellular Probes*, 16(3), 223–229. <https://doi.org/10.1006/mcpr.2002.0415>

Nie, K., Qi, S. xiang, Zhang, Y., Luo, L., Xie, Y., Yang, M. jie, Zhang, Y., Li, J., Shen, H., Li, Q., & Ma, X. jun. (2012). Evaluation of a Direct Reverse Transcription Loop-Mediated Isothermal Amplification Method without RNA Extraction for the Detection of Human Enterovirus 71 Subgenotype C4 in Nasopharyngeal Swab Specimens. *PLoS ONE*, 7(12). <https://doi.org/10.1371/journal.pone.0052486>

Notomi, T., Okayama, H., Masubuchi, H., Yonekawa, T., Watanabe, K., Amino, N., & Hase, T. (2000). Loop-mediated isothermal amplification of DNA. In *Nucleic Acids Research* (Vol. 28, Issue 12).

Yan, L., Zhou, J., Zheng, Y., Gamson, A. S., Roembke, B. T., Nakayama, S., & Sintim, H. O. (2014). Isothermal amplified detection of DNA and RNA. In *Molecular BioSystems* (Vol. 10, Issue 5, pp. 970–1003). Royal Society of Chemistry. <https://doi.org/10.1039/c3mb70304e>

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Francois, P., Tangomo, M., Hibbs, J., Bonetti, E. J., Boehme,

C. C., Notomi, T., Perkins, M. D., & Schrenzel, J. (2011). Robustness of a loop-mediated isothermal amplification reaction for diagnostic applications. *FEMS Immunology and Medical Microbiology*, 62(1), 41–48. <https://doi.org/10.1111/j.1574-695X.2011.00785.x>

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French industrial company specialized in particular in Life Sciences and Biodefense fields. Bertin Technologies is the Legal Manufacturer of the tests and markets these tests.

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Leading company in the field of veterinary infectious disease diagnostics. Enalees developed the RNA extraction technology used for the tests and the point-of-care test model.

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