



# AIRBORNE LATEX ALLERGENS CONTROL IN HOSPITAL ENVIRONMENTS

UNIPR (Italy)

## / CONTEXT

Latex products use have increased a lot over last 20 years with a considerable incidence on allergy to latex (single-use products). **Latex sensitivity in the healthcare environment** is an health problem, particularly, children poly-operated as Spina bifida patients and sanitary operators are concerned.

In this pilot study, we carried out the **environmental monitoring of latex allergen** in three different places in order to get quantitative data on airborne latex allergens in specific risky environments.

## / MATERIALS

- Coriolis® $\mu$ , sterile cones, 15 ml of sterile collection liquid.
- Traditional method on PTFE filters (25 mm  $\varnothing$  - 1 $\mu$ m).
- ELISA tests (Indoor Biotechnologies).

## / PROTOCOL

- Coriolis® $\mu$ : 2x20 min sampling for 3 days - 250 L/min - Quantification of liquid, sub-division, storage at -80°C.
- PTFE filters: 24 hours continuously sampling for three days -14 L/min - Allergen extraction overnight in PBS, centrifugation, sub-division, storage at -80°C.
- Specific ELISA tests for allergens: *Hev b1*, *Hev b3*, *Hev b5*, *Hev b6.02*.

## / CONCLUSION

These preliminary results indicate that Coriolis® $\mu$  is a suitable method for the sampling of airborne allergen of latex by reducing the time of sampling (20 min vs. 24h) and increasing the efficiency.

Furthermore, the Coriolis® $\mu$  air sampler collects any airborne particles and can also give information on specific airborne microorganisms, pollens, viruses...


## / RESULTS

- The results for *Hev b1* (ng/m<sup>3</sup>), *Hev b5* and *Hev b6.02* show the efficiency of the Coriolis air sampler to collect airborne latex allergens; *Hev b1* results are illustrated in the graph for the three rooms and for both equipments.

(1- a room of our laboratory in full operation and use of latex gloves, 2- an operating pediatric surgery room during normal activities, 2'- during activities in latex safe conditions and 3- near nursery room).

Coriolis vs. PTFE positive data :

- Room 1: 5+/6 vs. 2+/3
- Room 2: 2+/2 vs. 0+/1
- Room 2': 2+/2 vs. 0+/1
- Room 3: 2+/2



## / CUSTOMER



UNIVERSITÀ DEGLI STUDI DI PARMA