

Precellys 50 mL solution



Precellys 50 mL solution for large sample homogenization

- Exclusive vacuum-based 50 mL tube holder
- Easy and reliable solution for large sample homogenization
- Compatible with the Precellys Evolution

PRECELLYS 50 mL SOLUTION

Experience the reliability and efficiency of the Precellys Evolution for the 50 mL format specifically adapted to the needs of large sample homogenization. The dedicated 50 mL tube holder, developed for Precellys Evolution, ensures the efficient and uniform treatment of 1 to 3 x 50 mL tubes at the same time.

The Precellys 50mL tubes are compatible with the Precellys Evolution model and can be used with the dedicated 50mL tube holder.

Precellys 50 mL lysing kits (20 prefilled tubes):

50 mL Precellys Lysing kit VK01	P002475-LYSK0-A.0
50 mL Precellys Lysing kit CK28	P002478-LYSK0-A.0
50 mL Precellys Lysing kit CKMix50	P002476-LYSK0-A.0



50 mL tube holder pack for Precellys Evolution : S002474-PEVO0-A.0



Discover our comprehensive range of solutions



AIR
SAMPLING



SAMPLE
PREPARATION



ASSAY
KITS



CELL IMAGING
SYSTEM



TECHNOLOGY
ON DEMAND



Bertin Technologies SAS
Parc d'activités du Pas du Lac
10 bis, avenue Ampère
78180 Montigny-le-Bretonneux, FRANCE
01 39 30 60 00

WWW.BERTIN-TECHNOLOGIES.COM