

# Precellys 50mL solution



## Precellys 50mL solution for large sample homogenization

- Exclusive vacuum-based 50mL tube holder
- Easy and reliable solution for large sample homogenization
- Compatible with the Precellys® Evolution Touch\*

# PRECELLYS 50mL SOLUTION

Experience the reliability and efficiency of the Precellys® Evolution Touch for the 50mL format specifically adapted to the needs of large sample homogenization. The dedicated 50mL tube holder, developed for Precellys® Evolution Touch, ensures the efficient and uniform treatment of 1 to 3 x 50mL tubes at the same time.

The Precellys 50mL tubes are compatible with the Precellys® Evolution Touch\* model and can be used with the dedicated 50mL tube holder.

\* Also compatible with Precellys® Evolution

## Precellys 50mL lysing kits (20 prefilled tubes):

50mL Precellys Lysing kit VK01	P002475-LYSK0-A.0
50mL Precellys Lysing kit CK28	P002478-LYSK0-A.0
50mL Precellys Lysing kit CKMix50	P002476-LYSK0-A.0



## 50mL tube holder pack for Precellys Evolution Touch : S002474-PEVO0-A.0



Discover our comprehensive range of solutions



AIR  
SAMPLING

SAMPLE  
PREPARATION

ASSAY  
KITS

TECHNOLOGY  
ON-DEMAND



Bertin Technologies SAS  
Parc d'activités du Pas du Lac  
10 bis, avenue Ampère

78180 Montigny-le-Bretonneux, FRANCE  
01 39 30 60 00

[WWW.BERTIN-TECHNOLOGIES.COM](http://WWW.BERTIN-TECHNOLOGIES.COM)